

MicrUs and MicrUs Pro Series Ultrasound Systems

Echo Wave A Software

Measurements and Calculations Reference Manual



TELEMED
Ultrasound Medical Systems

Internet page: <http://www.pcultrasound.com>
Information E-mail: info@pcultrasound.com
Support E-mail: support@pcultrasound.com

Content

Echo Wave A Software Measurements and Calculations Overview.....	3
1 B mode general measurements and calculations.....	4
1.1 B Distance.....	4
1.2 B Area (Ellipse).....	4
1.3 B Volume (1 ellipse).....	5
2 M mode general measurements and calculations.....	6
2.1.1 Two-points M measurements.....	6
2.1.1.1 M Distance.....	6
2.1.1.2 M Time.....	7
2.1.1.3 M Velocity.....	7
2.1.1.4 M Heart Rate (HR).....	7
3 Revision History.....	8

Echo Wave A Software Measurements and Calculations Overview

This document presents equations that are used for Echo Wave A measurements and calculations.

- B mode measurements and calculations

Distance

Area (method: 1 ellipse)

Circumference(method: 1 ellipse)

Volume (method: 1 ellipse)

- M mode measurements and calculations

Distance

Time


Velocity

Heart Rate (method: 2 beats distance)

1 B mode general measurements and calculations

In this section are presented basic equations that are used both for general measurements and calculations. **Please note that not all here described measurements may have control items in software user interface, but they may be used in other calculations.**


1.1 B Distance



Distance d between points $pt1$ and $pt2$ is calculated using the following equation:

$$d(pt1, pt2) = \sqrt{(x1 - x2)^2 + (y1 - y2)^2} .$$

1.2 B Area (Ellipse)




During measurements and calculations we assume that ellipse axis between two circular marker points $pt1$ and $pt2$ is "long axis", and axis with one circular endpoint marker $pt3$ is "short axis". And this "long axis" and "short axis" notation remains unchanged no matter what are real lengths of these axes.

Area S and circumference P (perimeter) of an ellipse with long axis length $dl1$ and short axis length $ds1$ are calculated using the following equations:

$$S = \frac{\pi \cdot dl1 \cdot ds1}{4} , P = \pi \cdot \sqrt{\frac{1}{2}((dl1)^2 + (ds1)^2)} .$$

Volume V of an ellipsoid with axes lengths $d11$, $ds1$, and $ds2=ds1$ is calculated using the following equation:




$$V = \frac{\pi \cdot d11 \cdot ds1 \cdot ds2}{6}.$$

1.3 **B** Volume (1 ellipse)


See "B Area (Ellipse)" section.

2 M mode general measurements and calculations

Usually in M mode ultrasound image horizontal axis (x-axis) represents time (in seconds), and vertical axis (y-axis) represents depth (in millimeters). On M mode image are usually performed two-point -based measurements and calculations. For measurements and calculations we use (time [s], depth [mm]) coordinate system, where each point can be described by its time (in seconds [s]) and depth (in millimeters [mm]). For example, notation $pt1(x1,y1)=(5,120)$ means that coordinates of point $pt1$ are $x1=5s$ and $y1=120mm$.



2.1.1 Two-points M measurements



For two-point M measurements we use coordinates (time,depth) of two end-points $pt1$ and $pt2$ of one line (distance).

2.1.1.1 M Distance

Distance between points $pt1$ and $pt2$ is calculated using the following equation:

$$d = \text{abs} (d1 - d2),$$

here

- d [mm] - distance,
- $d1$ [mm] - depth at point $pt1$,
- $d2$ [mm] - depth at point $pt2$,
- $\text{abs}(\dots)$ means that is calculated absolute value.

2.1.1.2 *M Time*

Time interval (difference) between points pt1 and pt2 is calculated using the following equation:

$$t = \text{abs} (t1 - t2),$$

here

- t [s] - time interval (difference),
- t1 [s] - time at point pt1,
- t2 [s] - time at point pt2.

2.1.1.3 *M Velocity*

Velocity between points pt1 and pt2 is calculated using the following equation:

$$\text{Vel} = \text{abs} (d2 - d1) / \text{abs} (t2 - t1),$$

here

- Vel [mm/s] - velocity,
- t1 [s] - time at point pt1,
- d1 [mm] - depth at point pt1,
- t2 [s] - time at point pt2,
- d2 [mm] - depth at point pt2.

2.1.1.4 *M Heart Rate (HR)*

Heart Rate (HR) using markers pt1 and pt2 is calculated according to the following equation:

$$\text{HR} = 60 * \text{beats_num} / \text{abs}(t2-t1),$$

here

- HR [beats/min] - Heart Rate in beats per minute,
- abs(t2-t1) [s] - time interval between markers pt1 and pt2,
- beats_num [beats] - the number of heart beats (e.g., 2) in measured time interval.

3 Revision History

Revision	Revision Date	Description of Revision	Revision Author
1.0.0	2019.10.07	Initial Release	V.Perlibakas
1.0.1	2020.05.14	Changed first page photo.	V.Perlibakas
1.2.0	2020.08.11	Added information about M mode measurements.	V.Perlibakas
1.2.1	2020.08.13	Added section “M Heart Rate”.	V.Perlibakas